






Systemec D-200



- (1) Fuse ambient air ventilator
- (2) Fuse vacuum-pump
- (3) Fuse compressor
- (4) Fuse water-pump / valve
- (5) Cut-Off



- (6) Cooling water 3/4"
- (7) Demineralized water 3/4"
- (8) Compressed air
- (9) Interface RS 232
- (10) Power supply
- (11) Drain 3/4"



Legend

No.	Systemec DX	Systemec DE	Systemec DB	Remarks
1	X	X	X	
2	X	X	X	Vacuum-pump
3	X	X	X	
4	X	X	X	
5	X	X	X	
6	X	X	-	
7	X	X	-	
8	O	O	-	Only in conjunction with optional water cooling system
9	X	X	-	
10	380 - 400 V			Modified power supplies available on request
11	X	X	X	

X = Standard ; O = Optional erhältlich ; - = Not available

Allgemeine Technische Eigenschaften		Oberfläche		Maßstab	Seitenzahl
Größe	Material	Farbe	Textur	Verhältnis	von
Systemec				Systemec D-200en	1
Systemec				Systemec D-200en	1