






Systemec D-45



- (1) Condenser for compressor
- (2) Fuse ambient air ventilator
- (3) Fuse vacuum-pump
- (4) Fuse water-pump / valve
- (5) Fuse compressor
- (6) Cut-Off



- (7) Cooling water 3/4"
- (8) Demineralized water 3/4"
- (9) compressed air
- (10) Tank drain 1/4"
- (11) Drain 3/4"
- (12) Interface RS 232
- (13) Aqua-Stop
- (14) Power supply



Legend

No.	Systemec DX	Systemec DE	Systemec DB	Remarks
1	X	X	X	
2	X	X	X	
3	O	O	-	Only in conjunction with optional vacuum system
4	X	X	X	
5	X	X	X	
6	X	X	X	
7	X	X	-	
8	X	X	-	
9	O	O	-	Only in conjunction with optional water cooling system
10	X	X	X	
11	X	X	X	
12	X	X	-	
13	X	X	-	
14		220 - 230 V	-	Modified power supplies available on request

X = Standard ; O = Optional ; - = Not available

Allgemein: Referenz ISO 2908 11.04		Oberfläche Maßstab		Gewicht	
Systemec		Systemec D-45en		1	