




Systec V-65



A

- (1) Interface RS 232
- (2) Power supply

B

- (3) Fuse vacuum-pump
- (4) Fuse compressor
- (5) Fuse water-pump
- (6) Cut-Off

C

- (7) Compressed air
- (8) Demineralized water 3/4"
- (9) Cooling water 3/4"
- (10) Drain 3/4"

Legend

No.	Systec VX	Systec VE	Systec VB	Remarks
1	X	X	-	
2	380 - 400 V			Modified power supplies available on request
3	X	X	X	Vacuum pump
4	X	X	X	
5	X	X	-	
6	X	X	X	
7	O	O	-	Only in conjunction with optional water cooling system
8	X	X	-	
9	X	X	-	
10	X	X	X	

X = Standard ; O = Optional ; - not available

Allgemeine Information		Öberfläche	Maßstab	Gewicht
ISO 2566 11-m		Werkstoff:		
Datei		Reviz	Werkstoff	
Systec		Systec V65en		1
Systec		Systec V65en		AD