




Systemec V-75



- A**
- (1) Interface RS 232
 - (2) Power supply

- B**
- (3) Fuse vacuum-pump
 - (4) Fuse compressor
 - (5) Fuse water-pump
 - (6) Cut-Off

- C**
- (7) Compressed air
 - (8) Demineralized water 3/4"
 - (9) Cooling water 3/4"
 - (10) Drain 3/4"

Legend

No.	Systemec VX	Systemec VE	Systemec VB	Remarks
1	X	X	-	
2		380 - 400 V		Modified power supplies available on request
3	X	X	X	Vacuum pump
4	X	X	X	
5	X	X	-	
6	X	X	X	
7	O	O	-	Only in conjunction with optional water cooling system
8	X	X	-	
9	X	X	-	
10	X	X	X	

X = Standard ; O = Optional ; - not available

Allgemeine Technische Zeichnung	Oberfläche ISO 2768 TS max	Maßstab	Gewicht
Systemec V75en	Werkstoff:	1	A3