

	Scan[®] 100 Manual colony counter	Ref. 435 000
	Accurate view and work comfort. Scan [®] 100 is a high-tech manual colony counter remarkable for its user-friendly functions and wonderful ergonomomy. A LED lighting system with Dark-Field technology gives you an accurate and contrasted view of the colonies. The integrated USB port allows export of the results to guarantee traceability.	

TECHNICAL SPECIFICATIONS

- ▶ All stainless steel, very robust
 - ▶ Long lasting white LEDs, adjustable intensity
 - ▶ Dark Field lighting for a contrasted view of the colonies
 - ▶ Multifunction digital screen
 - ▶ For all Petri dishes: ø 55, 65, 90, 150 mm
 - ▶ Counting: from 0 to 1999 CFU
 - ▶ Reads on: PetriFilm[™], Sanita-kun[™]/RIDA[™]COUNT, Compact Dry[™] with multi-adapter (option)
 - ▶ Adjustable sensitivity of tactile surface
 - ▶ Adapted to any pen
 - ▶ Indicative sound: 4 levels
 - ▶ Ergonomic hand rest
 - ▶ Built-in accessories storage rack
 - ▶ Traceability, instant export of results via USB data transfer (Word, Excel[™], Notepad...)
-
- ▶ Power: ~100-240V / 50-60Hz
 - ▶ Machine: size (L x W x H): 25 x 30 x 18 cm, weight: 4.5 kg
 - ▶ Box: size (L x W x H): 50 x 40 x 40 cm, weight: 7.4 kg
-
- ▶ 1 year guarantee (after registration)
 - ▶ Life time service
 - ▶ Designed and made in France
- Manufactured under: ISO 9001 V2000, CE, WEEE 2002/96/EC, RoHS 2002/95/EC standards

DELIVERED WITH

- ▶ 1 white background plate
- ▶ ø 90 mm Petri dish adapter
- ▶ 1 black marker pen
- ▶ Power cable
- ▶ USB cable
- ▶ User manual, IQ/OQ/PQ guidelines
- ▶ Quick user guide

OPTIONS / ACCESSORIES

- ▶ Ref. 435 001 Flexible magnifying lens
- ▶ Ref. 435 002 Multi-adapter
- ▶ Ref. 435 010 Set of Spiral grids ø 90 and 150 mm (white and transparent grids)
- ▶ Ref. 435 020 Wolffhuegel grid